

SURGICAL BLADES & SCALPEL HANDLES



Product

Surgical blades are used for cutting skin and tissue during surgical procedures. The most common type of scalpel blade used in surgery is the non-attached, replaceable blade.

Features blades

- ✓ designed for tissue cutting and puncturing
- ✓ carbon steel material for exacting performance
- ✓ special edge treatment reduces drag and tissue trauma
- ✓ provide consistent sharpness – right out of the box
- ✓ fits all makes of handles - conforms to ISO 7740

Attributes handles

- ✓ slotted handle for better grip - provides greater strength and control
 - ✓ durable stainless steel material
 - ✓ swift lock fixation on handle
- weight 24 g
large 13,8 cm

Packaging

gamma sterilized
individually packaged in easy peel color coded foil for maximum safety
wrapped in aluminum V.C.I paper



**For quality
&
safety**




www.sutures.be

SURGICAL BLADES & SCALPEL HANDLES




For blade handle N° 3



Surgical blades without handle	ZB10	ZB11	ZB12	ZB12d	ZB15	ZB15c
Surgical blades with plastic handle	ZBD10	ZBD11	ZBD12		ZBD15	

For blade handle N° 4




Surgical blades without handle	ZB18	ZB20	ZB21	ZB22	ZB23	ZB24	ZB25
Surgical blades with plastic handle		ZBD20	ZBD21	ZBD22	ZBD23	ZBD24	

STITCH CUTTER



Unit	Order no
100 pc.	ZBS

SCALPEL HANDLES



Product	Unit	Order no
N° 3	10 pc.	ZBH3
N° 4	10 pc.	ZBH4

For quality
&
safety



www.sutures.be