

# ULTRA HIGH MOLECULAR WEIGHT POLYETHYLENE FIBER




## Product

E-FIBER is a new ultra-high tensile-strength suture for orthopaedic applications. Especially for treatment of ruptured Cranial Cruciate Ligament (CrCL) by extra-articular stabilisation technique.

## Features

- ✓ Unique braiding construction with a blue tracer
  - provides improved visibility
  - provides superior flexibility and tensile strength
- ✓ Offers durability, long lasting wear and fatigue resistance
- ✓ Non absorbable surgical suture
- ✓ Lubricious design allows each knot to slide for smoother
- ✓ More efficient performance and excellent knot security



For quality  
&  
safety

**smi**  
www.sutures.be

## E-FIBER SUTURES

		Unit	Order no.
metric 5 USP 2 length 500 mm	dogs < 10 kg	1 pc.	ZCB5
metric 6 USP 3&4 length 500 mm	dogs < 20 kg	1 pc.	ZCB6
metric 7 USP 5 length 500 mm	dogs < 30 kg	1 pc.	ZCB7
metric 9 USP 7 length 600 mm	dogs < 40 kg	1 pc.	ZCB9

## E-FIBER LATERAL SUTURES WITH STRAIGHT & CURVED NEEDLE

		Unit	Order no.
metric 5 USP 2 length 300 mm	dogs < 10 kg	1 pc.	ZCB5N
metric 6 USP 3&4 length 500 mm	dogs < 20 kg	1 pc.	ZCB6N
metric 7 USP 5 length 500 mm	dogs < 30 kg	1 pc.	ZCB7N
metric 9 USP 7 length 600 mm	dogs < 40 kg	1 pc.	ZCB9N

## E-FIBER LATERAL SUTURES WITH STRAIGHT & CURVED NEEDLE + 1 CRIMP

		Unit	Order no.
metric 5 USP 2 length 300 mm	dogs < 10 kg	1 pc.	ZCB5NC
metric 6 USP 3&4 length 500 mm	dogs < 20 kg	1 pc.	ZCB6NC
metric 7 USP 5 length 500 mm	dogs < 30 kg	1 pc.	ZCB7NC
metric 9 USP 7 length 600 mm	dogs < 40 kg	1 pc.	ZCB9NC

## E-FIBER LATERAL SUTURES DOUBLE LINE WITH STRAIGHT & CURVED NEEDLE

		Unit	Order no.
metric 9 USP 7 double line length 500 mm	dogs > 40 kg	1 pc.	ZCB9N2

## E-FIBER LATERAL SUTURES DOUBLE LINE WITH STRAIGHT & CURVED NEEDLE + 2 CRIMPS

		Unit	Order no.
metric 9 USP 7 double line length 500 mm	dogs > 40 kg	1 pc.	ZCB9N2C

