

ABSORBENT OP-COVER



Surgical drapes material

56 g/m² non breathable 3 layers
Polyolefines based film 24 g/m²
Rubber based adhesive (UFD technology) 2 g/m²
Hydrophilic polypropylene spunbond 30 g/m²

Colour

Blue

Attributes

- OP cover is latex free
- OP cover is soft and drape-able
- Liquid distributes and absorbs quickly

Product	Without hole	Order no
30 x 45 cm		ZOC3045
45 x 60 cm		ZOC4560
60 x 90 cm		ZOC6090
90 x 120 cm		ZOC90120
120 x 120 cm		ZOC120120
120 x 180 cm		ZOC120180
180 x 250 cm		ZOC180250

Product	Hole size	Order no
31 x 44 cm	1,5 x 4,5 cm	ZOC3144-S
31 x 44 cm	3,5 x 9 cm	ZOC3144-M
31 x 44 cm	5 x 14 cm	ZOC3144-B
44 x 74 cm	7 x 20 cm	ZOC4474-S
44 x 74 cm	8 x 24 cm	ZOC4474-M
44 x 74 cm	10 x 29 cm	ZOC4474-B



For quality
&
safety