

# SKIN STAPLER



## Product

Skin Staplers, also known as wound staplers or surgical staplers are medical tools for a variety of skin closures. The staples first penetrate the skin and are then formed, holding the tissue together.



## Benefits

- ✓ Excellent visibility of operating site
- ✓ Save considerable time for wound closure
- ✓ Non-stop feed action reduces jamming
- ✓ Better apposition of skin edges with much thinner scar after operation

## Features

- ✓ Contains 35 medical grade stainless steel staples (measure 6.4 mm x 4.0 mm)
- ✓ Pre-cock staple function for exact staple positioning
- ✓ Closed staple shaped squarely restricting movement or rolling of a staple
- ✓ Includes a counter to indicate the exact number of staples

### Staple size



Product	Unit	Order no
Skin stapler single, sterile packed	10	ZS35W

**For quality & safety**

# STAPLE REMOVER



## Product

SMI removers are specially designed and produced for the good function to remove the staples easily after the skin was sutured by our disposable skin stapler


## Features

- SMI removers will work with any brand of skin staple
- Fast and easy removal
- Reducing pain and infection risk

## Packing

Individually sterile packed

## How to remove



Product	Unit	Order no
Single use staple remover	20	ZSR
Single use staple remover with plastic handle	40	ZSR2

*For quality  
&  
safety*



[www.sutures.be](http://www.sutures.be)