


SUPRAMID

BLACK











Type	Twisted and coated multifilament Sizes 5/0 and 6/0 monofilament
Composition	Polyamide 6.6 strands enclosed in a polyamide 6 sheet
Colour	Black
Tissue reaction	Minimal
Absorption	Non absorbable, gradually encapsulated by connective tissue
Presentation	Box with 12 sutures







*For quality
&
safety*





smi




www.sutures.be


 Round bodied taper point		Metric	USP	Length	Order number
1/2 circle 	HR 16 mm 	1.5	4/0	75 cm	<i>black</i> 4150116
		2	3/0	75 cm	<i>black</i> 4200116
	HR 20 mm 	2	3/0	45 cm	<i>black</i> 4200120
		HR 26 mm 	3	2/0	75 cm


 Round bodied taper point		Metric	USP	Length	Order number
Straight 	GR 22 mm 	1.5	4/0	50 cm	<i>black</i> 4152022


 Reverse cutting		Metric	USP	Length	Order number	
1/2 circle 	HS 17 mm 	1.5	4/0	75 cm	<i>black</i> 4151417	
		2	3/0	75 cm	<i>black</i> 4201417	
	HS 20 mm 	1	5/0	75 cm	<i>black</i> 4101420	
		1.5	4/0	75 cm	<i>black</i> 4151420	
	HS 22 mm 	2	3/0	75 cm	<i>black</i> 4201420	
		1.5	4/0	75 cm	<i>black</i> 4151422	
	HS 26 mm 	2	3/0	75 cm	<i>black</i> 4201422	
		2	3/0	75 cm	<i>black</i> 4201426	
			3	2/0	75 cm	<i>black</i> 4301426

 Cutting		Metric	USP	Length	Order number
Straight 	GS 20 mm 	1.5	4/0	50 cm	<i>black</i> 4152120
	GS 25 mm 	1.5	4/0	60 cm	<i>black</i> 4152125

 Reverse cutting		Metric	USP	Length	Order number
3/8 circle 	DS 12 mm 	0.7	6/0	75 cm	<i>black</i> 4071512
		1	5/0	45 cm	<i>black</i> 4101512-45
		1	5/0	75 cm	<i>black</i> 4101512
		1.5	4/0	45 cm	<i>black</i> 4151512

Reverse cutting		Metric	USP	Length	Order number
 3/8 circle	DS 16 mm	0.7	6/0	75 cm	<i>black</i> 4071516
		1	5/0	45 cm	<i>black</i> 4101516-45
		1	5/0	75 cm	<i>black</i> 4101516
		1.5	4/0	45 cm	<i>black</i> 4151516
		1.5	4/0	75 cm	<i>black</i> 4151516-75
		2	3/0	45 cm	<i>black</i> 4201516
		2	3/0	75 cm	<i>black</i> 4201516-75
	DS 16 mm black	1	5/0	75 cm	<i>black</i> 4101516VB
	DS 19 mm	0.7	6/0	75 cm	<i>black</i> 4071519
		1	5/0	75 cm	<i>black</i> 4101519
		1.5	4/0	45 cm	<i>black</i> 4151519
		1.5	4/0	75 cm	<i>black</i> 4151519-75
		2	3/0	45 cm	<i>black</i> 4201519-45
		2	3/0	75 cm	<i>black</i> 4201519
	DS 19 mm Easy Pass	1.5	4/0	75 cm	<i>black</i> 4151519P-75
	DS 19 mm black	1	5/0	75 cm	<i>black</i> 4101519VB
		1.5	4/0	75 cm	<i>black</i> 4151519VB
	DS 24 mm	2	3/0	75 cm	<i>black</i> 4201524
		3	2/0	75 cm	<i>black</i> 4301524
		3.5	0	75 cm	<i>black</i> 4351524
	DS 36 mm	3	2/0	75 cm	<i>black</i> 4301536
3.5		0	75 cm	<i>black</i> 4351536	
4		1	75 cm	<i>black</i> 4401536	

Round bodied taper cutting point		Metric	USP	Length	Order number
 1/2 circle	HRT 18 mm	1.5	4/0	75 cm	<i>black</i> 4150518
		2	3/0	75 cm	<i>black</i> 4200518

Round bodied taper cutting point		Metric	USP	Length	Order number
 3/8 circle	DRT 18 mm	1.5	4/0	75 cm	<i>black</i> 4150418
		2	3/0	75 cm	<i>black</i> 4200418