

Material for repairing the cruciate ligaments

With this system you can unassisted and without much effort perform a solid and long lasting fixation of the cruciate repair.




CRUCIATE NEEDLE

				Order number
size 2	large	length 80 mm	dz.	ZC0000002
size 2	large	length 80 mm	6 pc.	ZC0000002-6
size 3	medium	length 62 mm	dz.	ZC0000003
size 3	medium	length 62 mm	6 pc.	ZC0000003-6
size 4	small	length 45 mm	dz.	ZC0000004
size 4	small	length 45 mm	6 pc.	ZC0000004-6
size 5	very small	length 30 mm	dz.	ZC0000005
size 5	very small	length 30 mm	6 pc.	ZC0000005-6



CRUCIATE REPAIR CRIMP TUBE

			Order number
10 mm (SS 316)		dz.	ZC00010
12 mm (SS 316)		dz.	ZC00012
14 mm (SS 316)		dz.	ZC00014



NYLON CRUCIATE SURGICAL SUTURES

			STERILE EO	Order number
CCL nylon sutures 50 lb x 500 mm	dogs < 15 kg	1 pc.		ZC050
CCL nylon sutures 80 lb x 800 mm	dogs < 20 kg	1 pc.		ZC080
CCL nylon sutures 100 lb x 800 mm	dogs > 20 kg	1 pc.		ZC100
	dogs > 40 Kg = 2 x 100 Lb			

NYLON CRUCIATE SURGICAL SUTURES WITH 1 CRIMP

			STERILE EO	Order number
CCL nylon sutures 50 lb x 500 mm + 10 mm crimp		1 pc.		ZC05010
CCL nylon sutures 80 lb x 800 mm + 12 mm crimp		1 pc.		ZC08012
CCL nylon sutures 100 lb x 800 mm + 12 mm crimp		1 pc.		ZC10012

NYLON CRUCIATE SURGICAL SUTURES WITH 3 CRIMPS

			STERILE EO	Order number
CCL nylon sutures 50 lb x 500 mm + 3 x 10 mm crimp		1 pc.		ZC05010-3
CCL nylon sutures 80 lb x 800 mm + 3 x 12 mm crimp		1 pc.		ZC08012-3
CCL nylon sutures 100 lb x 800 mm + 3 x 12 mm crimp		1 pc.		ZC10012-3

NYLON CRUCIATE SURGICAL SUTURES WITH NEEDLE + CRIMP

			STERILE EO	Order number
CCL nylon sutures 50 lb x 500 mm + very small 30 mm needle + 10 mm crimp		1 pc.		ZC0501005
CCL nylon sutures 80 lb x 800 mm + small 45 mm needle + 12 mm crimp		1 pc.		ZC0801204
CCL nylon sutures 100 lb x 800 mm + medium 62 mm needle + 12 mm crimp		1 pc.		ZC1001203

NYLON CRUCIATE SURGICAL SUTURES LOOP (DOUBLE LINE)

			STERILE EO	Order number
CCL nylon sutures 50 lb x 1000 mm + very small 30 mm needle loop		1 pc.		ZC0500005L
CCL nylon sutures 80 lb x 1600 mm + small 45 mm needle loop		1 pc.		ZC0800004L
CCL nylon sutures 80 lb x 1600 mm + medium 62 mm needle loop		1 pc.		ZC0800003L
CCL nylon sutures 100 lb x 1600 mm + medium 62 mm needle loop		1 pc.		ZC1000003L

NYLON CRUCIATE SURGICAL SUTURES WITH 2 CRIMPS LOOP (DOUBLE LINE)

			STERILE EO	Order number
CCL nylon sutures 50 lb x 1000 mm + very small 30 mm needle + 2 x 10 mm crimps loop		1 pc.		ZC0501005L
CCL nylon sutures 80 lb x 1600 mm + medium 62 mm needle + 2 x 12 mm crimps loop		1 pc.		ZC0801204L
CCL nylon sutures 100 lb x 1600 mm + medium 62 mm needle + 2 x 12 mm crimps loop		1 pc.		ZC1001203L